

Retort Clamps and Stands



1. BOSS HEAD

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body Double Mild Steel Thumb Type Screw Keys with Accepts Objects up to 16mm Dia Rod. Both Side

Cat No. / BH.001



2. BOSS HEAD

Made Of Pressure Die cast Aluminum Alloy, Powder coated Body with Double Mild steel Thumb Type screw key Accepts objects up to 16 mm Dia Rod. On one side & upto 19 mm Rod on the other side.

Cat No. / BH.002



3. BOSS HEAD

Made of Pressure Die cast zinc Alloy, Powder Coated Body with Double Mild steel Thumb Type screw Keys Accepts Objects up to 16 mm Dia Rod. Both Side

Cat No. / BH.003



4. BOSS HEAD

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body Double Mild Steel Thumb Type Screw Keys with Accepts Objects up to 16mm Dia Rod. Both Side

Cat No. / BH.004






5. BOSS HEAD RECTANGULAR

Made of Pressure Die cast Aluminum Alloy Powder Coated, Rectangular Shape. Double Mild steel Thumb Type Screw Keys Accepts Objects upto 18 mm Dia Rod. Both side.


Cat No. / BH.005



6. BOSS HEAD CAST IRON

Made of Cast iron Black Powder Coated Body with Double Mild steel Thumb Type screw Keys , Accepts Objects. Up to (A - 25 mm] (B -19 mm) (C - 15 mm] Dia Rod.


Cat No. / BH.006



7. BOSS HEAD CAST IRON

Made of Cast Iron Black Powder Coted Body With Brass Thumb Type Key Accepts Objects. Up to (A - 25 mm] (B -19 mm) (C - 15 mm] Dia Rod.

Cat No. / BH.007



8. BOSS HEAD BRASS


Mad of Brass Casting All Yellow Polished With Brass Thumb Type Screw Accepts Objects. Up to (A - 25 mm] (B -19 mm) (C - 15 mm] Dia Rod.

Cat No. / BH.008

9. BOSS HEAD BRASS

Made of Brass Casting All Nickel Platted With Brass Thumb Type Screw Accepts Objects. Up to (A - 25 mm) (B-19 mm) (C - 15 mm) Dia Rod.

Cat No. / BH.009



10. BOSS HEAD BRASS

Made Of Brass Casting Yellow Polished Body With Double Mild Steel Screw Key With Plastic Coated Tops, Accepts Objects, Upto 19mm Dia Rod.

Cat No. / BH.010



11. BOSS HEAD BRASS

Made Of Brass Casting Yellow Polished Body With Double Mild Steel Screw Key With Plastic Coated Tops, Accepts Objects, Upto 19mm Dia Rod.

Cat No. / BH.011



12. BOSS HEAD ROTATABLE

Made of pressure Die Cast zinc Alloy Powder coated, Rotatable Type, Lockable at any angle, Double Mild steel Thumb Type screw Keys or with plastic coated Tops. Accepts object up to 15 mm Dia Rod Both side.

Cat No. / BH.012





13. BOSS HEAD ROTATABLE

Made of Pressure Die cast zinc Alloy Powder coated, Rotatable Type. Lockable at any angle. Triple Mild steel screw Thumb Type Keys or Accepts object up to 15 mm Dia. on one side & 19 mm on the other side.

Cat No. / BH.013



14. BOSS HEAD ROTATABLE

Made of Pressure Die Cast Aluminum Alloy Powder Coated, Rotatable Type, Lockable at any angle, Double Mild Steel Screw Thumb Type Keys with or Plastic Coated Tops. Accepts Objects up to 19 mm Dia Rod. Both Side.

Cat No. / BH.014



15. RETORT CLAMP CROSS PATTERN EXTENTION

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CL.001



16. RETORT CLAMP CROSS PATTERN ROTATABLE

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Boss head with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CL.002



17. RETORT CLAMP CROSS PATTERN EXTENTION

Made of Pressure Die Cast Aluminium Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CL.003



18. RETORT CLAMP CROSS PATTERN ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CL.004





19. RETORT CLAMP CROSS PATTERN EXTENTION

Made of Sand Cast Aluminum Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminum Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CL.005



20. RETORT CLAMP CROSS PATTERN ROTATABLE

Made of Sand Cast Aluminum Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Rotatable Lockable at any angle. Fitted with a Bosshead with Mild Steel Screw Key Accepts Objects up to 100 mm Dia.


Cat No. / CL.006



21. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CL.007



22. THREE PRONG CLAMP ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Three Prong PVC Sleeves Covered Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CL.008

23. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Three Prong PVC Sleeves Covered Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.


Cat No. / CL.009



24. THREE PRONG CLAMP ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Three Prong PVC Sleeves Covered Jaw Clamp Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key. Accepts Objects up to 100 mm Dia.

Cat No. / CL.010



25. THREE PRONG CLAMP EXTENTION DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, Powder Coated, Three Prong PVC Sleeves Covered Jaw Clamp fitted to 6" or 9" Long Mild Steel Nickel Plated Rod. Jaw is Double Adjustable type with the help of Double Mild Steel Thumb Type Screw. Accepts Objects up to 65 mm Dia.

Cat No. / CL.011



26. THREE PRONG CLAMP ROTATABLE DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body. Three Prong PVC Sleeves Covered Jaw is Double Adjustable type with the help of Double Mild Steel Thumb Type Screw Keys. Clamp Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Keys. Accepts Objects up to 65 mm Dia.

Cat No. / CL.012





27. BURETTE CLAMP EXTENTION DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body Double "V" Shaped Cork Lined or PVC Covered Jaw Clamp fitted to 6" Or 9" Long Mild Steel Nickel Plated Rod. Jaw is Double Adjustable type with the help of Double Mild Steel Thumb type Screw Keys. Accepts Objects up to 65 mm Dia.

Cat No. / CL.013



28. BURETTE CLAMP ROTATABLE DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, Powder Coated, Double "V" Shaped Cork Lined or PVC Covered Jaw is Double Adjustable type with the help of Double Mild Steel Thumb Type Screw Keys. Clamp Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Screw Keys. Accepts Objects up to 65 mm Dia.


Cat No. / CL.014



29. BURETTE CLAMP EXTENTION

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body, Double "V" Shaped Cork Lined Jaw Clamp fitted to 6" Long Stainless Steel Pipe. Jaw is Adjusted with the help of fly nut which acts on spring action. Accepts Objects up to 60 mm Dia.

Cat No. / CL.015



30. BURETTE CLAMP ROTATABLE

Made of Pressure Die Cast Zinc Alloy, Powder Coated Body, Double "V" Shaped Cork Lined Jaw Clamp Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type screw Key. Jaw is Adjusted with the help of fly nut which acts on spring action. Accepts Objects up to 60mm Dia.

Cat No. / CL.016

31. THREE PRONG CLAMP EXTENTION DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, All Nickel Plated Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 6" or 9" Long Mild Steel Rod. Jaw is Double Adjustable type with the help of Double Mild Steel Screw Keys. Accepts Objects up to 65 mm Dia.

Cat No. / CL.017

32. THREE PRONG CLAMP ROTATABLE DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, All Nickel Plated Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted with Bosshead to 2" Long Mild Steel Rod, Clamp Rotatable & Lockable at any angle. Jaw is Double Adjustable type with the help of Double Mild Steel Screw Keys. Accepts Objects up to 65 mm Dia.

Cat No. / CL.018

33. BURETTE CLAMP EXTENTION DOUBLE KEYS


Made of Pressure Die Cast Zinc Alloy, All Nickel Plated Body, Double "V" Shaped Cork Lined or PVC Covered Jaw Clamp fitted to 6" or 9" Long Mild Steel Rod. Jaw is Double Adjustable type with the help of Double Mild Steel Screw Keys. Accepts Objects up to 65 mm Dia.

Cat No. / CL.019

34. BURETTE CLAMP ROTATABLE DOUBLE KEYS

Made of Pressure Die Cast Zinc Alloy, All Nickel Plated Body, Double "V" Shaped Cork Lined or PVC Covered Jaw is Double Adjustable type with the help of Double Mild Steel Screw Keys. Clamp fitted Bosshead to 2" Long Mild Steel Rod. Clamp Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Screw Key Accepts Objects up to 65 mm Dia.

Cat No. / CL.020





35. BURETTE CLAMP EXTENTION

Made of Pressure Die Cast Zinc Alloy, All Nickel Platted, Double "V" Shaped Cork Lined Jaw Clamp fitted to 6" Long Stainless Steel Pipe. Jaw is Adjusted with the help of fly nut which acts on spring action. Accepts Objects up to 60 mm Dia.

Cat No. / CL.021



36. BURETTE CLAMP ROTATABLE

Made of Pressure Die Cast Zinc Alloy, All Nickel Platted Body, Double "V" Shaped Cork Lined Jaw Clamp Rotatable & Lockable at any angle. Fitted with Bosshead mild steel screw Key. Jaw is Adjusted with the help of fly nut which acts on spring action. Accepts Objects up to 60 mm Dia.

Cat No. / CL.022



37. FOUR PRONG CLAMP EXTENTION

Made of Sand Cast Brass, Brass Polished Cork Lined four prong Jaw Clamp fitted to 9" & 12" long Mild Steel Chrome Plated Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.023



38. THREE PRONG CLAMP EXTENTION

Made of Sand Cast Brass, Brass Polished PVC Sleeve Covered Three Prong Jaw Clamp fitted to 9" & 12" long Mild Steel Chrome Plated Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.024

39. CONDENSER CLAMP EXTENTION

Made of Sand Cast Brass, Brass Polished Cork Lined Two Prong Jaw Clamp fitted to 9" & 12" long Mild Steel Chrome Plated Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.025

40. SUSPENSION CLAMP

Made of Pressure Die Cast Aluminum Alloy Powder Coated, Mild Steel "V" shape hook fitted in Aluminium Rod, Mild Steel Screw Key with Plastic Coated Top or Mild Steel Thumb Type Key. Useful for Suspending articles.

Cat No. / CL.026

41. RETORT CLAMP CROSS PATTERN EXTENTION


Made of Pressure Die Cast Zinc Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia. All Nickel Platted

Cat No. / CL.027

42. FOUR PRONG CLAMP EXTENTION

Made of Sand Cast Aluminum Alloy, Powder Coated Body, Four Prong Rubber Sleeves Covered Jaw Clamp fitted to 6" Inch Long Mild Steel Nickel Plated Rod. Accepts objects up 100 mm dia.

Cat No. / CL.028






43. FOUR PRONG CLAMP ROTATABLE

Made of Sand Cast Aluminum Alloy, Powder Coated Body, Four Prong Rubber Sleeves Covered Jaw Clamp fitted to Mild Steel Rod, Rotatable Lockable at any angle. Fitted with a Bosshead with Mild Steel Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CL.029



44. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted, Three Prong PVC Sleeves Covered Clamp fitted to 6" Inch Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts objects up 100 mm dia.

Cat No. / CL.030



45. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, PVC Sleeves Covered Jaw Clamp fitted to 6" Inch Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts objects up 100 mm dia.

Cat No. / CL.031



46. THREE PRONG CLAMP COPPER OXODIZED

Made of Sand Cast Brass, All Copper Oxodized, PVC Sleeve Covered Three Prong Jaw Clamp fitted to 9" long Mild Steel Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.032

47. FOUR PRONG CLAMP COPPER OXODIZED

Made of Sand Cast Brass, All Copper Oxodized, Cork Lined Jaw Clamp fitted to 9" long Mild Steel Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.033



48. CONDENSER CLAMP COPPER OXODIZED

Made of Sand Cast Brass, All Copper Oxodized, Cork Lined Jaw Clamp fitted to 9" long Mild Steel Rod. Accepts objects up to 100 mm Dia.

Cat No. / CL.034



49. RETORT CLAMP CROSS PATTERN EXTENTION

Made of Pressure Die Cast Alumimium Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CL.035



50. RETORT CLAMP CROSS PATTERN ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Four Prong Cork Lined Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CL.036





51. TRIPOD STAND

Cast Iron Circular Top with concentric circles fitted on Mild Steel Zinc Plated Rods. Available Sizes: Ht. x Top Dia.- (8"x 6")

Cat No. / TPS.001



52. TRIPOD STAND

Cast Iron Circular Top fitted on Mild Steel Zinc Plated Rods. Available Sizes:Ht. x Top Dia.-(A - 6x4) (B - 7x4) (C - 7x5) (D - 8x5)

Cat No. / TPS.002



53. TRIPOD STAND

Cast Iron Triangular Top fitted on Mild Steel Zinc Plated Rods. Available Sizes: Ht. x Top Dia.- (A - 6x4) (B - 7x4) (C - 7x5) (D - 8x5)

Cat No. / TPS.003



54. TRIPOD STAND STAINLESS STEEL

Stainless Steel Wire Triangular Top fitted on Stainless Steel Rods. Available Sizes: Ht. x Top Dia.- (6"x 4") (8"x 5")


Cat No. / TPS.004

55. TRIPOD STAND STAINLESS STEEL

Stainless Steel Wire Circular Top fitted on Stainless Steel Rods.

Available Sizes: Ht. x Top Dia.- (6"x 4") (8" x 5")

Cat No. / TPS.005



56. TRIPOD STAND

Cast Iron Circular Top fitted on Mild Steel Zinc Plated Rods.

Available Sizes: Ht. x Top Dia.- (A - 3.5"x6") (B - 4"x6") (C - 5"x8") (D - 6"x8")

Cat No. / TPS.006





57. RETORT RING EXTENTION

Mild Steel Nickel Plated Extension Ring.
Available Sizes: 50mm, 75mm, 100mm, 125mm, 150mm

Cat No. / RR.001



58. RETORT RING WITH BOSS HEAD

Mild Steel Wire Retort Ring fitted with Bosshead.
Available Powder Coated & Nickel Platted Sizes:
50mm, 75mm, 100mm, 125mm, 150mm

Cat No. / RR.002



59. RETORT RING

Cast Iron Retort Ring with Bosshead.
Available Sizes: 50mm, 75mm, 100mm, 125mm, 150mm

Cat No. / RR.003



60. RAY BOX

Made of Pressure Die Cast Aluminum Alloy,
Powder Coated Body With Two Slide Mild Steel
Sheet With Two Lenses With Holder.

Cat No. / RB.001

61. CRUCIBLE TONG

Made of Stainless Steel Wire Fine Finishing. In All Size Available.


Cat No. / CT.001



62. CRUCIBLE TONG

Made of Stainless Steel Wire Fine Finishing. In All Size Available.

Cat No. / CT.002



63. BEAKER TONG

Made of Stainless Steel Wire Fine Finishing. In All Size Available.

Cat No. / BT.003



64. FLASK TONG

Made of Stainless Steel Wire Fine Finishing. In All Size Available.

Cat No. / FT.004





65. TEST TUBE HOLDER

Made of mild Steel & Brass wire Nickel Plated.

Cat No. / TTH.001



66. SLIDE TRAY

Made of Aluminium Sheet Anodised.


Cat No. / ST.001



67. PENDULUM PULLEY

Made of Aluminium Streep with ABS Moulded Round Pulley, Singal, Double & Triple Pulley .

Cat No. / PP.001



68. TEST TUBE STAND

Made of Aluminium sheet all sizes Available.

Cat No. / TTS.001

69. RETORT STAND

Retort Stand Base made from mild steel sheet. Bent from 2 sides, Center hole on the smaller side. Powder coated body. Available Sizes: 8" x 5", 9" x 6", 12" x 6"

Cat No. / RS.001



70. RETORT STAND

Retort Stand Base made from mild steel sheet. Bent from all sides Center hole on the smaller side. Powder coated heavy base.

Available Sizes : (6" x 4") (8" x 5") (9" x 6") (10" x 6") (11" x 6")

Cat No. / RS.002



71. RETORT STAND

Retort Stand Base made from mild steel sheet. Bent from 2 sides, Center hole on the smaller side. Powder coated body.

Available Sizes: (6" x 4") (7" x 5") (8" x 5") (9" x 6") (10" x 6") (12" x 6")

Cat No. / RS.003




72. RETORT STAND

Retort Stand Base made from mild steel sheet. Bent from all sides, Center hole on the smaller side. Powder coated heavy base with rubber feet on two sides for superior gripping.

Available Sizes : (8" x 5")

Cat No. / RS.004





73. RETORT STAND

Retort Stand Base made from cast iron. Center hole on the smaller side.

Available Sizes: (6" x 4") (7" x 5") (8" x 5") (9" x 6") (11" x 6") (12" x 6") (12" x 7")

Cat No. / RS.005



74. GAS TAP

Mad of Brass Casting Fine Finishing all Chrome Plated Air Regulator With Stop Cock Available, One way, Two way, Three way, Four way,

Cat No. / GT.001



75. RETORT ROD ALUMINIUM ANODISED

Anodised Aluminum Rod for Retort Stands. Available All Sizes : As per requirement.

Cat No. / RRA.001



76. RETORT ROD MILD STEEL CHROME PLATED

Chrome plated & Zinc Plated Mild Steel Rod for Retort Stands. Available All Sizes: As per requirement.

Cat No. / RRMS.002

77. RETORT ROD STAINLESS STEEL

Stainless Steel Rod for Retort Stands. Available All Sizes : As per requirement.

Cat No. / RRSS.003



78. BURETTE CLAMP DOUBLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Powder Coated, Spring Loaded Jaws for (2) Double Burette.

Cat No. / BC.001



79. BURETTE CLAMP SINGLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Powder Coated, Spring Loaded Jaws for (1) Single Burette.

Cat No. / BC.002



80. BURETTE CLAMP DOUBLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Nickel Plated Body, Spring Loaded Jaws for (2) Double Burette.

Cat No. / BC.003





81. BURETTE CLAMP SINGLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Nickal Plated Body, Spring Loaded Jaws for (1) Single Burette.

Cat No. / BC.004



82. BURETTE CLAMP DOUBLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Spring Loaded Jaws for (2) Double Burette.

Cat No. / BC.005



83. BURETTE CLAMP SINGLE

Burette Clamp made of Pressure Die Cast Aluminum Alloy, Powder Coated Body, Spring Loaded Jaws for (1) Single Burette.

Cat No. / BC.006



84. COMPASS

Made of Aluminium Sheet Body One Side.

Cat No. / CP.001



85. COMPASS TRANSPARENT

Made of Aluminium Sheet Body Both Side.

Cat No. / CP.002



86. SPRIT LAMP STAINLESS STEEL

Made of StainLess Steel Sheet.

Cat No. / SL.001



87. SPRIT LAMP GLASS

Made of Glass.

Cat No. / SL.002





88. B.P. BULB

Cat No. / BPB.001



89. ARMLET BAG

Cat No. / AB.001



90. LABORATORY TUBING

Cat No. / LT.001



91. CORRACE

Cat No. / CE.001

92. BUNSEN BURNER WITH STOP COCK

Made of Cast Iron Powder Coated Base With Brass Pipe Chrome Platted, Screw Type Air Regulator With Brass Stop Cock.

Cat No. / BB.001

93. BUNSEN BURNER WITH NOZZEL

Made of Cast Iron Powder Coated Base With Brass Pipe Chrome Platted With Nozzle.

Cat No. / BB.002

94. BUNSEN BURNER DOUBLE FLEME

Made of Cast iron Powder Coated Base With Brass Pipe Chrome Platted, Screw Type Air Regulator With Brass Stop Cock.

Cat No. / BB.003

95. BUNSEN BURNER DOUBLE FLEME

Made of Cast iron Powder Coated Base With Brass Pipe Chrome Platted, With Nozzle.

Cat No. / BB.004





96. BOSS HEAD

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted body Double Mild Steel Thumb Type Screw Keys with Accepts Objects up to 16mm Dia Rod. Both Side

Cat No. / BHP.001



97. BOSS HEAD

Made Of Pressure Die cast Aluminum Alloy, Nickel platted body with Double Mild steel Thumb Type screw key Accepts objects up to 16 mm Dia Rod. On one side & up to 19 mm Rod on the other side.

Cat No. / BHP.002



98. BOSS HEAD

Made of Pressure Die Cast Zinc Alloy, Nickel Platted Body Double Mild Steel Thumb Type Screw Keys with Accepts Objects up to 16mm Dia Rod. Both Side.

Cat No. / BHP.003



99. BOSS HEAD


Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body Double Mild Steel Thumb Type Screw Keys with Accepts Objects up to 16mm Dia Rod. Both Side.

Cat No. / BHP.004

100. BOSS HEAD RECTANGULAR

Made of Pressure Die cast Aluminum Alloy Nickel Platted, Rectangular Shape. Double Mild steel Thumb Type Screw Keys Accepts Objects up to 18 mm Dia Rod. Bothside.

Cat No. / BHNP.005



101. BOSS HEAD ROTATABLE

Made of Pressure Die cast zinc Alloy Nickel platted, Rotatable Type. Lockable at any angle. Triple Mild steel screw Thumb Type Keys or Accepts object up to 15 mm Dia. on one side & 19 mm on the other side.


Cat No. / BHNP.006



102. BOSS HEAD ROTATABLE

Made of Pressure Die Cast Aluminum Alloy Nickel platted, Rotatable Type, Lockable at any angle, Double Mild Steel Screw Thumb Type Keys. Accepts Objects up to 19 mm Dia Rod. Both Side.

Cat No. / BHNP.007



**103. RETORT CLAMP CROSS PATTERN EXTENTION**

Made of Pressure Die Cast Zinc Alloy, Nickel Platted Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.001**104 . RETORT CLAMP CROSS PATTERN ROTATABLE**

Made of Pressure Die Cast Zinc Alloy, Nickel Platted Body, Four Prong Cork Lined Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Boss head with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.002**105. RETORT CLAMP CROSS PATTERN EXTENSION**

Made of Pressure Die Cast Alumimium Alloy, Nickel Platted Body, Four Prong Cork Lined Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.003**106. RETORT CLAMP CROSS PATTERN ROTATABLE**

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Four Prong Cork Lined Jaw Clamp fitted to 2" Long Mild Steel Nicel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.004

107. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.005



108. THREE PRONG CLAMP ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.006



109. THREE PRONG CLAMP EXTENTION

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.007



110. THREE PRONG CLAMP ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Three Prong PVC Sleeve Covered Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.008





111. RETORT CLAMP CROSS PATTERN EXTENSION

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Four Prong Cork Line Jaw Clamp fitted to 6" Long (Mild Steel Nickel Plated / Aluminium Anodised) Rod. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.009



112. RETORT CLAMP CROSS PATTERN ROTATABLE

Made of Pressure Die Cast Aluminum Alloy, Nickel Platted Body, Four Prong Cork Line Jaw Clamp fitted to 2" Long Mild Steel Nickel Plated Rod, Rotatable & Lockable at any angle. Fitted with a Bosshead with Mild Steel Thumb Type Screw Key. Accepts Objects up to 100 mm Dia.

Cat No. / CLNP.010

