

Tapered Cork Stopper



MEDILAB offers the highest quality of regular, short, extra long, extra extra long tapered corks in agglomerated quality.



Sizes Available

We provide tapered cork stoppers in any size and per customer's requirement



Tapered Cork sizes and measurements						
Stopper #	Lenght		Top Dia.		Top Dia.	
	Inch	MM	Inch	MM	Inch	MM
0000	1/2	12.7	3/16	4.8	3/32	2.4
000	1/2	12.7	1/4	6.3	5/32	4.0
00	1/2	12.7	5/16	7.9	7/32	5.5
0	1/2	12.7	3/8	9.5	9/32	7.1
1	5/8	15.9	7/16	11.1	21/64	8.3
2	11/16	17.5	1/2	12.7	3/8	9.5
3	3/4	19.0	9/16	14.3	27/64	10.7
4	13/16	20.6	5/8	15.9	15/32	11.9
5	7/8	22.2	11/16	17.5	17/32	13.5
6	15/16	23.8	3/4	19.0	37/64	14.7
7	1	25.4	13/16	20.6	5/8	15.9
8	1-1/16	27.0	7/8	22.2	43/64	17.1
9	1-1/8	28.6	15/16	23.8	47/64	18.6
10	1-1/4	31.8	1	25.4	49/64	19.4
11	1-1/4	31.8	1-1/16	27.0	53/64	21.0
12	1-1/4	31.8	1-1/8	28.6	57/64	22.6
13	1-1/4	31.8	1-3/16	30.2	61/64	24.2
14	1-1/4	31.8	1-1/4	31.8	1-1/64	25.8
15	1-1/4	31.8	1-5/16	33.3	1-5/64	27.4
16	1-1/2	31.8	1-3/8	34.9	1-3/32	27.9
17	1-1/2	31.8	1-7/16	36.5	1-5/32	29.4
18	1-1/2	31.8	1-1/2	38.1	1-7/32	30.9
19	1-1/2	31.8	1-9/16	39.7	1-9/32	32.5
20	1-1/2	31.8	1-5/8	41.3	1-11/32	34.1
21	1-1/2	31.8	1-11/16	42.9	1-13/32	35.7
22	1-1/2	31.8	1-3/4	44.5	1-15/32	35.7
23	1-1/2	31.8	1-13/16	46.0	1-17/32	38.9
24	1-1/2	31.8	1-7/8	47.6	1-19/32	40.5

Natural Tapered Cork Stopper

Wide at the top and narrow at the bottom, naturally tapered cork stoppers are extracted from the bark of *Quercus suber*. The fact that it's natural makes it a premium option for usage in the wine, beverages, pharmaceutical, homeware, aromatics and spice industry.



Sizes Available

We provide natural tapered cork stopper in all international sizes from size no. 0000 to size no. 24.

Agglomerated Tapered Cork Stopper

The agglomerated tapered cork stoppers or composition cork stoppers are usually produced from cork granules bound to each other by synthetic agents. This is where it differs from the naturally tapered cork stoppers.



Sizes Available

We provide agglomerated tapered cork stopper in all international sizes from size no. 0000 to size no. 24.

Customized Size & Bung Tapered Cork Stopper

Besides its usage in all other industries, bungs have found recognition in the pottery and glass jar segment. This is primarily not because the natural aspects of the corks compliment that of the ceramic, clay and glass works. But because the tapered corks fit the unequal openings of ceramic and glass containers with absolute precision.

Corks are not just natural but have a unique characteristic of moulding itself to the width of the variable openings of ceramic or clay containers. This is the reason behind cork's triumph over plastic or any other rigid material in this industry segment.

Besides bung tapered cork stoppers, bark-top cork stoppers are also used largely these days for both functional and decorative purposes. With the unprocessed natural chunk of bark of the mother tree left on the top of the tapered bung, it is usually used as a lid in containers used to store dry substances.

The bark-top and Bung cork stoppers are available in various sizes.



Customised Shape Cork Stopper

As the name suggests, these cork stoppers are custom-made to suit the usage and requirements of the end-users. By using advanced technology, 100% natural cork is given the due shape and size to churn out a piece exactly as requested by the customers.

Corks are increasingly replacing the usage of plastic and other artificial materials in the wine, beverages, pharmaceutical, homeware, sports, aromatics and spice industries. With the advantage of customisation coming into place, it makes far easier for these industries to obtain cork stoppers of their requirements hasslefree.

- ◆ Laboratories
- ◆ Chemical industry
- ◆ Pharmaceutical firms / homeopathy
- ◆ Olive oil
- ◆ Vinegar
- ◆ Decoration
- ◆ Cork floats
- ◆ Glassware
- ◆ Pottery



Cork Stand

Cork rings or cork stands, produced from multiple grades of cork material, are used in laboratories to hold glass flasks. At Advance Cork these cork stands are available in different sizes and shapes as per customers' requirements. These cork stands also come handy to support glass decorative flasks or vases at homes.

Size	Fits Flasks (ml)	Ld. (mm)	o.d. (mm)
1	10-100	30	80
2	200-500	60	110
3	1,000- 3,000	90	140
4	5,000	120	170
5	12,000	150	210



Features:

- ◆ Flawless concentric bores to easily support the round bottom of glass flasks
- ◆ Perfectly crafted parallel sides for aesthetic feel
- ◆ Beveled-finish on the inside to let the flasks sit comfortably
- ◆ Cork stands come with 30mm height to avoid the flasks to touch the ground
- ◆ Defined standards are adequately followed to produce these stands using quality cork

