

RUBBER LABWARE

SILICON RUBBER TUBING					
Catalog No.	Inner Dia ± 0.5mm	Wall ± 0.5mm	Outer Dia ± 0.5mm		
19808/01	1	1	3		
19808/02	2	1	4		
19808/03	3	1	5		
19808/04	3	1.5	6		
19808/05	4	1	6		
19808/06	4	1.5	7		
19808/07	4	2	8		
19808/08	5	1	7		
19808/09	5	1.5	8		
19808/10	6	1	8		
19808/11	6	1.5	9		
19808/12	6	2	10		
19808/13	6	3	12		
19808/14	7	1.5	10		
19808/15	7	2	11		
19808/16	7	3.5	14		
19808/17	8	1.5	11		
19808/18	8	2	12		
19808/19	8	3	14		
19808/20	9	2	13		
19808/21	10	2	14		
19808/22	10	3	16		
19808/23	12	2	16		
19808/24	12	3.5	19		
19808/25	14	3	20		
19808/26	15	4	23		
19808/27	16	3	22		
19808/28	25	3	31		
*Available in Roll of 10 Meters					



RUBBER LABWARE

SILICON RUBBER PLATINUM CURED TUBING

Catalog No.	Inner Dia ± 0.5mm	Wall ± 0.5mm	Outer Dia ± 0.5mm
19809/01	3	1.5	6
19809/02	4	1.5	7
19809/03	5	1.5	8
19809/04	6	1.5	9
19809/05	6	2	10
19809/06	7	1.5	10
19809/07	7	2	11
19809/08	8	1.5	11
19809/09	8	2	12
19809/10	9	2	13
19809/11	10	2	14
19809/12	12	2	16

*Available in Roll of 10 Meters



PRESSURE TUBING GREEN

Catalog No.	Inner Dia ± 0.5mm	Wall ± 0.5mm	Outer Dia ± 0.5mm
19806/01	6.4	3.2	12.8

BLOOD PRESSURE TUBE GREEN SUPERIOR

Catalog No.	Inner Dia ± 0.5mm	Wall ± 0.5mm	Outer Dia ± 0.5mm
19807/01	3	2.5	8



RUBBER LABWARE

LABORATORY RUBBER TUBING					
Catalog No.	Bore ± 0.5mm	Wall ± 0.5mm			
19803/01	4	1			
19803/02	5	1.5			
19803/03	6	1.5			
19803/04	7	1.5			
19803/05	8	1.5			
19803/06	9	2			
19803/07	10	2			
19803/08	12	2			
LABORATORY PRESSURE RUBBER TUBING					
Catalog No.	Bore ± 0.5mm	Wall ± 0.5mm			
19804/01	3	3			
19804/02	5	3			
19804/03	6	3			
19804/04	8	3			
19804/05	9	3			
19804/06	10	4			
19804/07	12	4			
LABORATORY PRESSURE HIGH RUBBER TUBING					
Catalog No.	Bore ± 0.5mm	Wall ± 0.5mm			
19805/01	3	5			
19805/02	5	5			
19805/03	6	5			
19805/04	8	5			
19805/05	9	5			
19805/06	10	6			
19805/07	12	6			







RUBBER LABWARE					
LABORATORY RUBBER CORKS					
Catalog No.	No.	Bottom ± 0.5mm	Top ± 0.5mm	Height ± 0.5mm	
19813/01	000	6	9	16	
19813/02	00	8	10.5	20	
19813/03	0	9	13	16	
19813/04	1	11	14	18.5	
19813/05	2	13	16	20	
19813/06	3	14	19	21	
19813/07	4	15	21	22	
19813/08	5	17.5	22.5	23	
19813/09	6	20	24.5	24	
19813/10	7	20	26	27.5	
19813/11	8	21	26.5	27.5	
19813/12	9	21.5	28	31	
19813/13	10	26	30	33	
19813/14	11	29	33	32	
19813/15	12	30	34.5	34.5	
19813/16	13	31	38.5	37	
19813/17	14	32	41.5	42.5	
19813/18	15	38	42	40.5	
19813/19	16	35	45	36.5	
19813/20	17	40	49.5	40.5	
19813/21	18	48	53.5	33.5	
19813/22	19	44	58	45	
19813/23	20	51	60	56.5	
19813/24	21	47	62	40	
19813/25	22	51	65	40	
19813/26	23	68	74	42	






RUBBER LABWARE					
SILICONE RUBBER CORKS					
Catalog No.	No.	Bottom ± 0.5mm	Top ± 0.5mm	Height ± 0.5mm	
19816/01	5	5	7	16	
19816/02	6	6	8	16	
19816/03	8	8	10.5	20	
19816/04	9	9	11.5	20	
19816/05	11	11	14	24	
19816/06	13	13	16	24	
19816/07	15	15	28	24	
19816/08	17	17	20	26	
19816/09	18	18	21	26	
19816/10	19	19	22.5	28	
19816/11	21	21	24.5	28	
19816/12	23	23	25.5	28	
19816/13	25	25	28.5	28	
19816/14	27	27	31	32	
19816/15	29	29	33	32	
19816/16	31	31	36	35	
19816/17	33	33	38.5	38.5	
19816/18	35	35	40.5	38.5	
19816/19	37	37	42.5	38.5	
19816/20	39	39	44.5	38.5	
19816/21	41	41	47	42	
19816/22	43	43	49	42	
19816/23	45	45	51	42	
19816/24	49	49	55	52	






RUBBER LABWARE					
RUBBER CONE FOR FILTRATION FLASK					
Catalog No.	Size				
19814/01	250/500ml				
19814/02	1000/2000ml				
19814/03	5ltr				
19814/04	10ltr				
RUBBER CONE FOR FILTRATION FLASK					
Catalog No.	Size				
19814/05	250/500ml				
19814/06	1000/2000ml				
19814/07	5ltr				
19814/08	10ltr				
RUBBER CONE (GUKO) (SET OF 7Pc)					
Catalog No.	Rubber Type				
19819/01	Natural Rubber				
19819/02	Silicone Rubber				
RUBBER SEPTA WITH TURN OVER FLANGE(PLAIN/SERRATED)					
Catalog No.	Size				
19822/01	12mm				
19822/02	15mm				
19822/03	18mm				
19822/04	25mm				
19822/05	B14				
19822/06	B19				
19822/07	B24				
19822/08	B29				





RUBBER LABWARE					
RUBBER SEPTA WITH TURN OVER FLANGE(PLAIN/SERRATED)					
Catalog No.	Size				
19821/01	12mm				
19821/02	15mm				
19821/03	18mm				
19821/04	25mm				
19821/05	B14				
19821/06	B19				
19821/07	B24				
19821/08	B29				
RUBBER BULBS					
Catalog No.	SIZE (Approx Capacity)				
19823/01	SMALL(5ml)				
19823/02	MEDIUM(10ml)				
19823/03	LARGE(15ml)				
19823/04	EXTRA LARGE(20ml)				
19823/05	OVAL(25ml)				
19823/06	SUPER LARGE(50ml)				
19823/07	EXTRA SUPER LARGE(100ml)				
RUBBER BULBS					
Catalog No.	SIZE (Approx Capacity)				
19823/08	SMALL(5ml)				
19823/09	MEDIUM(10ml)				
19823/10	LARGE(15ml)				
19823/11	EXTRA LARGE(20ml)				
19823/12	OVAL(25ml)				
19823/13	SUPER LARGE(50ml)				
19823/14	EXTRA SUPER LARGE(100ml)				









RUBBER LABWARE						
RUBBER TEATS						
Catalog No.	SIZE(Approx Capacity)	Rubber Type				
19826/01	1ml(CONICAL)	Natural Rubber				
19826/02	2ml(CONICAL)	Natural Rubber				
19826/03	1ml(CONICAL)	Silicone Rubber				
19826/04	2ml(CONICAL)	Silicone Rubber				
19826/05	SMALL	Latex Rubber				
19826/06	BIG	Latex Rubber				
AUTOMIZER BULB						
Catalog No.	Size					
19825/01	Small					
19825/02	Big					
PIPETTE FILLER						
Catalog No.	Rubber Type					
19827/01	Naturall Rubber					
19827/02	Silicone Rubber					
SILICONE CORK FILTRATION APPARATUS						
Catalog No.	Size					
19836/01	Small					
19836/02	Big					

RUBBER LABWARE						
VACCUM TUBE						
Catalog No.	Size					
19837/01	6mm					
19837/02	8mm					
LAB RUBBER SHEET						
Catalog No.	Rubber Type					
19834/01	Naturall Rubber 9.5"x9.5"x2.7mm					
19834/02	Silicone Rubber 9.5"x9.5"x2.7mm					
LAB RUBBER MAT						
Catalog No.	Rubber Type					
19832/01	Naturall Rubber 12"X12"					
19833/01	Silicone Rubber 12"X12"					
FEET						
Catalog No.	Size					
19838/01	Small					
19838/02	Medium					
19838/03	Large					

RUBBER LABWARE						
CRUCIBLE CONES						
Catalog No.	Size					
19839/01	30ML					
19839/02	50ML					
SILICONES RIDER FOR VOLUMETRIC FLASKS						
Catalog No.	SIZE	For Flask Capacity				
19840/01	1	5ml to 25ml				
19840/02	2	50ml				
19840/03	3	100ml to 250ml				
19840/04	4	500ml				
19840/05	5	1ltr				
RUBBER BASE FOR ROUND BOTTOM FLASK						
Catalog No.	Capacity					
19841/01	25ml					
19841/02	50ml					
19841/03	100ml					
19841/04	250ml					
19841/05	500ml					
19841/06	1ltr					
19841/07	2ltr					
RUBBER BASE FOR ROUND BOTTOM FLASK						
Catalog No.	Capacity					
19842/01	25ml					
19842/02	50ml					
19842/03	100ml					
19842/04	250ml					
19842/05	500ml					
19842/06	1ltr					
19842/07	2ltr					

RUBBER LABWARE						
SILICONE WEIGHT RINGS						
Catalog No.	Size					
19843/01	SMALL					
19843/02	BIG					
LATEX BLADDER FOR GAS SAMPLING						
Catalog No.	Size					
19844/01	NO.5					
19844/02	DOUBLE NOZZLE					
Other Rubber Labware Items						
Catalog No.	Item Name					
19828	Bunsen Burner Reinforced Tube (1 mtr. Orange)					
19830	Silicone Hand Protector					
19831	Silicone Flask Stand					
19845	Silicone Policeman					

RUBBER LABWARE						
19846	Silicone Spreader					
19847	Silicone Washer For Burette					
19848	Rubber Hammer					
19849	Rubber Ball For Tuning Fork					
19850	Rubber Pad					
19851	Rubber Tuning Fork Activator					
19852	Corrugated Sheet 300x150cm					
19853	Corrugated Tube 102cm					